

Using Custom Icon

RELEASE 2.6.0 OR ABOVE

You can set customized icons for Synap Editor.

Related API

- [static: addIcons](#)
- [static: getIcons](#)

Setting customized icons

- Add customized icons using the `SynapEditor.addIcons(Object)` API.
- Set the icons in the form of an `` or `<svg>` tag on an object whose function name is a key.
 - You can check for a detailed list of names of the available functions in the toolbar settings.
 - (It should be set before creating Synap Editor as new.)

```
<script>
    var customIcons = {
        'bold': '',
        'italic': '<svg xmlns="http://www.w3.org/2000/svg" version="1.1" viewBox="0 0 315 315"
>< polygon style="fill:#005ECE;" points="280,40 280,0 140,0 140,40 180.955,40 91.228,275 35,275 35,315 175,315
175,275 134.045,275 223.772,40 "></svg>'
    };
    SynapEditor.addIcons(customIcons); // Setting customized icons

    window.editor = new SynapEditor('synapEditor', synapEditorConfig); // Create Synap Editor
after setting.
</script>
```

